

Assut Medical Sàrl

At the service of surgery

Swiss - Quality

Assut Medical | product chart



Absorbable sterile sutures











AssuCryl®

Polyglycolic acid (PGA)

Braided, coated

Polycaprolactone & Calcium stearate

8-0 to 3&4 violet 6-0 to 2 undyed

60-90 days

65% at 14 days 40% at 21 days

Excellent handling
Excellent knot security
Very low capillarity

5 years

1 dozen or 3 dozen boxes

AssuCryl® Lactin

Polyglactin (PGLA)

Braided, coated

Poly(glycolide-co-L-Lactide) + calcium stearate (Glacomer 37 or 30/70)

> 8-0 to 3&4 violet 6-0 to 2 undyed

> > 56-70 days

75% at 14 days 50% at 21 days

Excellent handling
Excellent knot security
Very low capillarity

5 years

1 dozen or 3 dozen boxes

AssuCryl® Rapid

Polyglycolic acid (PGA)

Braided, coated

Polycaprolactone & Calcium stearate

6-0 to 2 undyed

40-45 days

50% at 6-9 days 0% at 17-21 days

Excellent handling
Excellent knot security
Very low capillarity

5 years

1 dozen or 3 dozen boxes

AssuCryl® MonoSlow

Polydioxanone (PDO)

Monofilament

None

7-0 to 2 violet

180-210 days

80% at 14 days 50% at 28 days

High initial tensile strength Excellent pliability Smooth tissue passage

5 years

1 dozen or 3 dozen boxes

AssuCryl® MonoRapid

Poliglecaprone (PGCL)

Monofilament

None

6-0 to 2 violet 6-0 to 2 undyed

90-120 days

50% at 6-8 days 0% at 21 days

Excellent pliability
Smooth tissue passage
Excellent knot security

5 years

1 dozen or 3 dozen boxes

At the service of surgery

Assut Medical | product chart



Non-absorbable sterile sutures











Astralen®	Silk	Nylon	Supramid	Polypropylene
Polyester	Silk	Polyamide	Polyamide	Polypropylene
Braided, coated	Braided, coated	Monofilament	Pseudo-monofilament or monofilament	Monofilament
Silicone	Wax, silicone & natural gum	None	Polyamide 6 for pseudo-monofilament None for monofilament	None
7-0 to 7 green or white	7-0 to 5 black 8-0 to 7-0 black virgin 9-0 to 7-0 blue virgin	7-0 to 3&4 blue 10-0 to 3&4 black	6-0 to 6 black	10-0 to 2 blue 7-0 to 2 clear
Smooth tissue passage Excellent knot security Easy to handle	Smooth tissue passage Excellent knot security Easy to handle	Effortless passage through tissue Excellent tensile strength	Superior flexibility Excellent knot security Smooth tissue passage	Effortless tissue passage Flexible and ties securely Avoids capillarity
5 years	5 years	5 years	5 years	5 years
1 dozen, 1/2 dozen or 3 dozen boxes	1 dozen or 3 dozen boxes	1 dozen or 3 dozen boxes	1 dozen or 3 dozen boxes	1 dozen or 3 dozen boxes

At the service of surgery

Assut Medical | product chart



Non-absorbable sterile sutures





AssuTopFiber®

Ultra high molecular weight Polyethylene & Polyester

Co-braided

Uncoated

2-0 to 7 green or white 2-0 to 7 white/blue or white/green

Outstanding tensile strength
Extremely strong
Flexible and smooth

5 years

1 dozen or 3 dozen boxes

Surgical Steel

Stainless steel

Monofilament

None

3&4 to 7 natural

Very flexible High tensile strength and unlikely to break

5 years

1 dozen long pack boxes

Specialty Products





Astralen® Tape

Polyester

Braided

None

3 or 5 or 7mm green or white

Easy to handle
Excellent knot security
Extremely strong

5 years

1 dozen or 1/2 dozen boxes

PTFE Pledgets

PTFE felt

Compressed fibers

None

Pre-cut in round, oval or rectangular, undyed

Pre-punched Used in combination with suture material

5 years

1 dozen or 3 dozen boxes

At the service of surgery

Assut Medical | contact





Business office

Assut Medical Sàrl Avenue de Lavaux 35 1009 Pully Switzerland

T +41 21 728 64 35 F +41 21 728 64 49



Factory address

Assut Medical Sàrl Succursale Sur le Crêt 13 2606 Corgémont Switzerland

T +41 32 489 31 70 **F** +41 32 489 31 76

@ info@assutsutures.com





Headquarters

Assut Medical Sàrl P.O Box 5 Av. de Rochettaz 57 1009 Pully Switzerland

www.assutsutures.com